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**Instruction Manual
Homecare
Patient Lift
AGILE**

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
Please notice:

In the text of these instructions square brackets and bold figures = [00] refer to the parts of the patient lift.

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Burmeier GmbH & Co. KG

PREFACE

Dear customer,
BURMEIER would like to thank you for the confidence you have placed in us and our products in deciding to purchase this AGILE homecare patient lift.

Each appliance has been tested by the manufacturer for electrical safety and functionality and has left our factory in perfect condition.

This instruction manual informs you as the operator and your users about all the functions necessary to ensure ease of operation and safe handling of this appliance on a daily basis.

You should, therefore, also regard this instruction manual as a practical reference book to be kept near the patient lift and at hand at all times.

We wish you and your users every success in operating the appliance. We are confident that this product enables us to play an important role in achieving just that.

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GENERAL NOTES

Before putting the patient lift into service for the first time:

■ Read through this instruction manual from start to finish so as to prevent any damage occurring due to incorrect operation.

Before using a patient lift, the user must check that the appliance is fully functional and in perfect working order and must follow the instruction manual in accordance with the Medizinprodukte-Betreiberverordnung (German abbreviation: MPBetreibV, Operators of Medical Products Ordinance) § 2. The same applies for accessories.

This patient lift fulfils all the requirements of the 93/42/EEC Medical Products Directive. It is classified as a Class 1 active medical product in accordance with the Medizinproduktegesetz (German abbreviation: MPG § 13, Medical Products Act).

The patient lift has been type-tested by an independent testing institute. As with any electrical device, incorrect use can prove hazardous.

Comply with your obligations as the owner under the Medical Product Use Directive to ensure that this medical product can be operated safely at all times without being a danger to patients, users or third parties.

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This instruction manual contains safety information which must be followed. All users working on and with the patient lift must be familiar with the contents of this instruction manual and follow the safety information provided.

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INSTRUCTIONS FOR THE OPERATOR

■ Any piece of technical equipment, electrical or otherwise, can prove hazardous if not properly operated. You are obligated to instruct users in the proper use of this patient lift in accordance with MPBetreibV § 5 (Operators of Medical Products Ordinance).

■ Ensure that users know where this instruction manual is located, in accordance with MPBetreibV § 9 (Operators of Medical Products Ordinance).

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DEFINITIONS OF INVOLVED PERSONS

In this manual the following persons are involved:

Operator

An operator is any natural or legal person who uses the patient lift or on whose instruction it is used (e.g. medical suppliers, specialized retailers, health insurance companies). It is the responsibility of the operator to instruct the user appropriately.

User

Users are persons who as a result of their vocational training, experience or briefing are authorized to operate the patient lift or to carry out work on it, or are instructed in handling the lift. Furthermore the user can recognise and avoid potential dangers and assess the clinical condition of the patient.

Patient

In this instruction manual, a patient is defined as a person who is infirm, in need of care or disabled and who is transferred with the patient lift.

Specialist staff

Qualified personnel are employees of the operator who as a result of their vocational training or briefing are entitled to deliver and to make the patient lift ready for operation.

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SAFETY INSTRUCTIONS

At the time of leaving the factory, this patient lift represents state-of-the-art technology and has been tested by an independent testing institute.

■ The most important objective of the safety information is to prevent personal injuries.

MEANING OF THE SAFETY SYMBOLS

In these instructions the following safety symbols are used:

Warning about injuries to persons

Dangerous voltage. Life threatening.

General danger. Injury or life hazard.

Warning about property damages

Possibility of damage to drive, material or environment.

Other symbols

Useful tip. For easier operation or better understanding of the unit.

The safety symbol used does not replace the text of the safety note. Therefore read the safety instructions and follow them exactly!

SAFETY INSTRUCTIONS FOR OPERATORS

■ Using this instruction manual, which must be provided with the patient lift, ensure that every user is instructed in the safe operation of this appliance before using it for the first time. Draw every user's attention to the possible hazards that can arise if the patient lift is improperly used. This applies in particular to the use of electrical actuators.

■ Patient lifts are classified as Class 1 active medical products according to § 13 of the German Medical Products Act (MPG). In order to ensure the permanently safe operation of this medical product, with no risk to patients, users or third parties, you must observe your obligations in accordance with the Operators of Medical Products Ordinance (MPBetreibV)! If the patient lift is in long-term use, test the functions and check for visible damage (see Chapter 6) after an appropriate period of time (recommendation: annually).

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- Only permit persons who have been properly instructed to use this patient lift.
- Make sure that substitute staff are also sufficiently well instructed in the safe operation of the patient lift!

- Ensure that the lifting accessories (harnesses, straps, slings etc.) are securely attached or hooked in.

If lifting accessories are used (seat slings, padded harnesses or straps), care must always be taken to securely hook in the devices used.

In use, the swivel bar attachment mainly lies transversely to the axis of the lift arm so that the patient can be lifted up or set down.

- When the patient is lifted up, ensure that fabric parts or straps are not caught in the wheelchair or bed.
- If the patient is turned in the lift, this operation is not performed by the swivel bar attachment but only by the person being lifted. Care must be taken to ensure that patients' limbs do not collide with metal parts as a result of any strong patient movements (e.g. by spastic persons).
- Since the swivel bar attachment of the patient lift must accommodate lateral play as well a full 360° rotation, care must always be taken when patients are turned.

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If you have any questions or concerns, please contact BURMEIER.

- Check that your staff is complying with the safety information.

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SAFETY INSTRUCTIONS FOR USERS

- Ensure that the operator instructs you in the safe operation of this patient lift.
- Each time before using the patient lift, check that it is in perfect working order.
- In order to avoid damage, ensure that no obstacles could impede adjustments to the patient lift.
- The "Agile" patient lift is designed for use in buildings (mainly in the homecare sector).

Since only temperatures of above 20°C and below 35°C occur in this type of location, the possibility of persons being exposed to high-temperature metal parts (max. 41°C) can be ruled out.

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It must be ensured that the lift is not permanently exposed to sunlight, in order that the surface temperature and the ambient temperature of the electrical components are maintained within the admissible range.

It is assumed that the controls and the control components are not exposed to hazards due to external force fields or electromagnetic emissions.

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You will find a checklist for evaluating the proper condition of the patient lift in Chapter 6 of this instruction manual.

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■ When not in use, always stow the handset in such a way that it cannot fall to the floor (hang it on the hook).

■ Make sure that the cable cannot be damaged by moving parts of the patient lift.

■ The adjustments may then only be carried out by, or in the presence of, a person instructed in the proper operation of the patient lift.

■ Ensure that seat slings and harnesses and straps are hooked in securely.

■ Make sure that seat slings and harnesses and straps are used correctly.

■ Check that the seat sling or harness or strap used is correctly positioned.

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PRODUCT DESCRIPTION

NORMAL OPERATING PROCEDURE

■ The "AGILE" patient lift is a lifting device for handicapped persons. This lift is mainly for use for transfers in the homecare sector (e.g. wheelchair, bed, toilet etc.).

Padded harnesses, straps and seat slings etc are used as lifting accessories.

■ The "AGILE" is designed as a traditional circular arc lift. The pedal-operated spreadable legs allow the patient lift to fit round wide wheelchairs.

■ The low profile construction of the legs also allows access underneath low beds.

■ Due to the double-tubed construction of the lift mast, the lift arm has an extremely high side load stability.

■ The swivel bar attachment of the patient lift is available in various sizes. The unit is controlled using a hand control with a coiled cable.

■ The power is provided by a rechargeable battery pack.

■ The control box with an integrated battery is fitted with an electronic overload protection and has a soft start/stop function.

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■ The control box has two light-emitting diodes (LEDs)

- ◆ Yellow LED: external battery charger is connected
- ◆ Green LED: handset is activated

■ An EMERGENCY OFF switch, and an electrical and a mechanical emergency lowering device, guarantee safe operation.

■ A sticker is provided on the lift arm giving the admissible safe working load.

This patient lift may only be operated under the operating instructions laid out in this instruction manual.

Any other use is deemed to be inappropriate.

DESCRIPTION OF MATERIALS

For the most part, the patient lift is manufactured from steel profiles whose surfaces are finished with a powder coating.

All surfaces are recognised as being safe for contact with skin.

STRUCTURAL DESIGN

Base and lift mast

- The pedal-operated spreadable legs are fitted with a foot support.
- The front castors are \varnothing 75 mm dual castors to allow for easy low-level access.
- The locking rear castors are 125 mm in diameter, which is ideal for handling the patient lift.
- The lift mast contains:
 - ◆ The lift arm actuator
 - ◆ Two push handles
 - ◆ The complete electrical power supply.

Electrical adjustment system

- The electrical adjustment system of this patient lift is fault-secure and flame-retardant (VO) and is made up of:
 - ◆ A control box with integrated EMERGENCY OFF and electrical emergency lowering device as well as an integrated battery.
 - ◆ A battery charger.
 - ◆ The electric motors for the lift actuator.
 - ◆ A handset with a stable hook.

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ASSEMBLY AND START UP PROCEDURE

■ The patient lift is supplied in three sections.

- 1) Base
- 2) Lift mast
- 3) Swivel bar attachment

■ The lift is assembled at the place of use through the qualified staff of the lift operator.

■ In view of the easy handling characteristics only one person is required to assemble the lift.

ASSEMBLY

☞ Remove the base from the packaging.

☞ Lock both rear castors.

☞ Take out the patient lift mast and insert this into the round tube socket in the base of the lift.

☞ Align the lift mast and secure the unit in place by inserting the two hexagon socket screws using an Allen key.

☞ Insert the round bolts for the swivel bar attachment into the tube socket for the lift arm and secure the unit in place with spring bolts.

INITIAL OPERATION

■ We recommend charging the battery for approx. 18 hours before using the lift for the first time to compensate for any loss of capacity during storage.

■ During the charging process the EMERGENCY OFF switch must **not** be activated.

■ During the charging process it is not possible to operate the lift electrically.

■ The battery must be charged at an ambient temperature of 0°C to 40°C. The most effective charging temperature range is between 5°C and 35°C. The capacity available is influenced by the temperature and the discharge current.

■ Ensure that the EMERGENCY OFF switch is activated. (By rotating the switch slightly the red emergency off button pops out).

☞ Using the handset, move the actuator in the UP direction and back again.

■ The patient lift is now ready for operation.

■ For the charging procedure, the battery charger supplied must be inserted into the bottom left of the control box.

The actuator system is fitted with an acoustic alarm. If a button is pressed the acoustic "charge signal" (recharge battery) sounds.

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SPECIAL SAFETY INSTRUCTIONS FOR THE DRIVE UNIT

■ If the load is too high, an electronic overload switch is activated and the control unit is automatically switched off.

Once the overload has been removed, the actuator system is reactivated when buttons are pressed again.

■ Continuous operation must not exceed **5 minutes**.

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HAND SWITCH


The electrical lift functions can be activated by the patient or the user on the handset.

■ The direction of travel of the actuators is shown by arrows / pictograms

■ The electric motors operate as long as the corresponding buttons are pressed.

■ The handset can be hung at any position on the patient lift with an elastic hook.

■ The coiled cable provides ample freedom of motion.

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CASTORS


The patient lift is fitted with two \varnothing 75 mm dual castors at the front and two \varnothing 125 mm rear castors. Both rear castors have brakes.

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EMERGENCY LOWERING OF THE LIFT ARM

■ In case of failure of the control box or a disconnection in the handset cable the lift arm can be lowered by inserting a ballpoint pen or similar object into the Emergency hole in the control box.

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■ In a completely current-free state the lift arm can be lowered manually using the quick release facility which is located at the upper end of the actuator rod.

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CLEANING AND DISINFECTING

GENERAL

Cleaning is the most important requirement for a successful chemical disinfection.

Generally a routine cleaning of the lift, if used by the same patient, is hygienically sufficient.

Disinfection of the lift is only necessary in case of visible contamination with infectious or potentially infectious material (blood, stool, pus) or by presence of an infectious disease under the direction of a physician.

The lift must be cleaned and washed down with disinfectant in case of a patient change!

Before beginning the cleaning, please note:

- Make sure that all plugs are properly inserted.
- There must be no external damage visible on any of the electrical components. Non-compliance with this advice could lead to malfunctions or damage to the electrical components as a result of penetrating water or cleaning solutions.
- The electrical components must not be cleaned with a water jet, a high pressure cleaner or other similar devices! Clean only with a moist cloth!
- If you suspect that water or any other form of moisture has penetrated the electrical components, take the patient lift out of service. Report this to the operator immediately.
- Failure to follow this safety advice could result in considerable damage to the equipment and lead to subsequent malfunctions!

CLEANING AND DISINFECTION PLAN

- Remove any straps and attached parts which do not belong to the lift.
- Clean all surfaces with a mild and environmentally friendly cleaning agent. Clean the handset in the same manner.
- If the patient lift has been visibly contaminated with infectious or potentially infectious materials, the patient lift should be subsequently disinfected by wiping with one of the disinfection media approved by the DGHM (Deutsche Gesellschaft für Hygiene und Mikrobiologie, German Society for Hygiene and Microbiology) which is suitable for the corresponding surfaces. The same applies for all lifts with patients who have notifiable diseases according to § 6 of the Infektionsschutzgesetz (IfSG, Protection against Infection Act), bacterial infections, or infections

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with multiple-resistant pathogens (e.g. MRSA, VRE), as well as all lifts in intensive care stations and infectious disease clinics. For all disinfections, the concentrations given in the DGHM list must be observed.

■ Disinfection of the castors is only necessary if they have been visibly contaminated with infectious or potentially infectious materials.

Note: continuous disinfection is only necessary in hospitals when a patient has a multiple-resistant pathogen (e.g. MRSA).

INSTRUCTING THE USER AND QUALIFIED PERSONNEL

In order to ensure that cleaning and disinfection are properly conducted, we recommend that users and experts are appropriately instructed. When providing instruction, observe the following points:

■ The clean lift must be transported to the patient's home by a suitable means to prevent it from becoming dirty or contaminated.

■ When the lift is dismantled, it must be cleaned and disinfected by wiping without delay. Staff should be informed of the special measures required for cleaning and disinfection and should carry out the procedure in a reliable manner (the operator should specify the operating procedure or the individual procedural steps). Care must be taken that only disinfection agents approved by the DGHM (German Society for Hygiene and Microbiology) are used, and that these are used only in the DGHM approved concentrations. The disinfection agent must be suitable for use with the surfaces to be disinfected.

■ For this activity, the staff should be outfitted with disposable aprons and gloves which are impermeable to fluids.

■ For the cleaning treatment, only fresh, clean cloths may be used which are subsequently sent to the laundry service.

■ After the cleaning treatment, and before leaving the cleaning location, the staff should disinfect their hands before proceeding to the next location. The staff should be equipped with a suitable pump dispenser containing a disinfection medium for hands.

By cleaning the lift directly in the place of use, it is possible to avoid clean lifts coming into contact with "unclean" lifts when transported in the same van. In this way, the transfer of potentially infectious germs, which may be found on the used patient lift, is prevented. A transfer of germs in terms of a nosocomial infection can be safely avoided by consistently and thoroughly following these recommendations.

If the lift is not to be immediately reused, it should be stored (covered) such that it is protected from dust, inadvertent dirt and contamination.

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CLEANING SUPPLIES AND DISINFECTANTS

Pay attention to the following recommendations to ensure that the functions and usability of the patient lift are preserved as long as possible:

- The surfaces must be intact; each sign of damage must be promptly repaired.
- We recommend (moist) cleaning by wiping. When selecting cleaning agents, ensure that the ones chosen are mild (gentle on skin and surfaces) and environmentally friendly. A standard household cleaner can generally be used.
- For the cleaning and care of coated metal parts, a moist cloth with a mild, commercial household cleaner is best.
- For disinfection by wiping, the disinfection medium used should be one which is listed by the DGHM (at the recommended concentration) and approved for the intended cleaning. Media and concentrations as given in the list from the Robert Koch Institute (or the former Bundesgesundheitsamt, Federal Health Office) must only be used in the case of an epidemic, in accordance with the Infektionsschutzgesetz (Protection against Infection Act), and only by order of the Health Office responsible.
- Despite the excellent mechanical resistance of the patient lift, scratches, markings, etc., which permeate the entire coating should be resealed using a suitable medium to prevent the penetration of moisture. For further information, consult BURMEIER or a specialist of your choice.

- Do not use scouring agents, stainless steel cleaners, abrasive cleaning agents or scouring pads. These substances can damage the surfaces.

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Note:

As a rule, aldehyde-based disinfection media have the advantage that they have a wide range of impact, a relatively low protein effect and are environmentally friendly. The main disadvantage of these agents is their potential to cause allergies and irritation. Glucoprotamine-based formulations do not have this disadvantage and are equally effective although most are somewhat more expensive. Disinfection media based on compounds which could potentially release chlorine may be corrosive for metals, synthetics, rubbers and other materials over longer contact periods or when concentrations are too high. Furthermore, these media have a higher so-called protein effect, are mucous membrane irritants and demonstrate poor environmental compatibility.

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CONTACT WITH DISINFECTANTS

- Pay attention to the exact dosage! We recommend the use of automated dosing instruments.
- Always prepare solutions with cold water in order to avoid the formation of vapours which are mucous membrane irritants.
- Wear gloves, in order to avoid direct skin contact.
- Do not keep prepared surface disinfection solutions in open containers with floating cleaning cloths. Be sure to cover all containers!
- Use sealable bottles with pump dispensers for moistening the cleaning cloths.
- Ventilate the room after the disinfection has been completed.
- Disinfect by wiping; do not disinfect by spraying! When spraying, a large portion of the disinfection medium is released as spray and could be inhaled.
- Furthermore, the wiping effect plays a significant role.
- Do not use alcohol for the disinfection of large surfaces.

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Legal Principles

In accordance with:

- ◆ Medizinprodukte-Betreiberverordnung § 4 (German Operator's Ordinance on Medical Products, Maintenance)
- ◆ Berufsgenossenschafts-Vorschrift BGV A3 (Directive of the German Employers Liability Insurance Association), previously BGV A2 (Testing of mobile electrical equipment in industrial use)

the operator is obligated to preserve the safe operating condition of medical products throughout their entire service life. This also includes professional servicing, to be carried out at regular intervals, and regular safety checks.

Recommendation

All electrical and mechanical components must be inspected once a year. In addition, the charging cable and the handset cable must be inspected each time that they have been subjected to mechanical strain and after each move to a new location.

The manufacturer assumes no liability whatsoever if the authorised service company or health supplies shop does not use original spare parts or accessory components.

The checklists on the following pages can be used as an aid for checking purposes:

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Checklist: Test protocol by the user

Test		OK	Not OK	Defect description
Visual check of electrical components				
Stickers and type plates	Present, legible			
Control unit housing, Battery pack	Correct position, damage			
Housing and actuator tube of motor	Correct position, damage			
Handset	Damage, membrane			
Motor cable, handset cable, power supply cable, battery pack	Damage, routing			
Plug of control unit housing	Available, correct position			
Visual check of mechanical components				
Stickers and type plates	Present, legible			
Lift mast	Damage, deformations			
Disconnect mast from lift base	Check stop angle			
Lift arm	Damage, deformations			
Lift arm	Split welded seams			
Swivel bar attachment	Damage			
Castors	Fixed position, damage			
Base legs	Damage, screw connection tight			
Base	Split welded seams			
Wearing parts such as joints	Damage			
Function test of electrical components				
End-of-travel cut-out of the motors	Automatic cut-out			
Handset	Function test			
Battery charger	Function test			
Control unit and motors	Noise level			
Control unit and motors	Function test			
For rechargeable battery systems:				
Electrical emergency lowering	Function test			
Locking device	Function test			
Integrated battery charger	Function test			
Signature of the tester:	Result of the test:			Date:

- If damage or a malfunction is suspected, the lift must be withdrawn from service immediately until the defective parts have been repaired!
- Please contact the operator responsible if a replacement or repair is necessary.

The test protocol draft on **pages 20 and 21** can be used.

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Customer / med. Facility / practice:			
Address:			
Carried out: <input type="checkbox"/> Repeat inspection		<input type="checkbox"/> Inspection prior initial operation	
<input type="checkbox"/> _____		<input type="checkbox"/> Inspection following repairs / servicing	
Equipment type: <input checked="" type="checkbox"/> patient lift	<input type="checkbox"/> Lifting aid	Protection classification: <input type="checkbox"/> I	<input checked="" type="checkbox"/> II
Lift type: <i>AGILE Patient Lift</i>	Inventory Number:		
Location:	Serial Number:		
Manufacturer: <i>Burmeier GmbH & Co. KG</i>	User-specific: <i>none</i>		
Testing equipment used (Type / Inventory Number):	1.		
MPG Classification: <i>Class I</i>	2.		
I. Visual Inspection	OK	Not OK	Defect Description
Visual inspection of the electrical components			
Stickers and type plates	Present, legible		
Control unit housing, battery pack	Correct position, damage		
Housing and actuator tube of motor	Correct position, damage		
Handset	Damage, membrane		
Motor cable, handset cable, mains cable, battery pack	Damage, routing		
Plug and plug cover of control unit	Available, correct position		
Visual inspection of the mechanical components			
Stickers and type plate	Present, legible		
Lift mast	Damage, Deformations		
Disconnect mast from base	Check stop angle		
Lift arm	Damage, Deformation		
Lift arm	Split welded seams		
Swivel bar attachment	Damage		
Castors	Damage, Fixed position		
Base legs	Damage, Screw connection tight		
Base	Split welded seams		
Wearing parts such as joints	Damage		

III. Performance check		OK	Not OK	Defect Description
Performance check of the electrical components				
End-of-travel cut-out of the motors	Automatic cut-out			
Handsets	Function test			
Battery charger	Function test			
Control unit and motors	Noise level			
Control unit and motors	Function test			
For rechargeable battery systems:				
Electrical emergency lowering	Function test			
Rechargeable battery locking facility	Function test			
integrated battery charger	Function test			
Performance check of the mechanical components				
Joints and pivots	Smooth operation			
Lift arm, emergency lowering	Function test			
Swivel bar attachment	Function test			
Castors	Braking, running			
Lift mast locking device	Function test			
Leg spreader, rocker pedal	Function test			
Leg spreader with max. load	Function test			
Lift mast with max. load	Function test			
Accessories (e.g. harnesses, straps, scales)	Fastening, damage, suitability			
Result of the inspection:				
All values within permissible range: <input type="checkbox"/> yes <input type="checkbox"/> no		Next inspection date:		
Inspection passed: <input type="checkbox"/> yes <input type="checkbox"/> no				
If inspection was not passed: <ul style="list-style-type: none"> <input type="checkbox"/> Defective, do not use the lift! => Repair <input type="checkbox"/> Defective, do not use the lift! => Take out of Service <input type="checkbox"/> Lift does not meet the safety standards 				
Test approval sticker applied: <input type="checkbox"/> yes <input type="checkbox"/> no				
Part of this inspection report:				
Comments:				
Inspected on:	Inspector:	Signature:		
Evaluated on:	Operator / Expert:	Signature:		

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SPARE PARTS

The corresponding replacement parts can be obtained from BURMEIER by specifying the item number, the order number and the serial number. The necessary details are found on the type plate which is located at the head end of the patient lift frame.

In order to maintain operational safety and the right to claim under warranty, only original BURMEIER replacement parts may be used!

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MANUFACTURER'S ADDRESS

For ordering replacement parts, customer service requests and other questions, please contact:

Burmeier GmbH & Co. KG
(A subsidiary of the Stieglmeyer Group)

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REPLACEMENT OF ELECTRICAL COMPONENTS

Dangerous voltage, electric shock possible

- Before doing any work unplug the main power cable from the wall socket!
- Any work on the electrical equipment may only be performed by customer service, the motor manufacturer or qualified and authorized electricians under consideration of the applicable VDE rules and safety guidelines!

- All the drive components are maintenance-free and must not be opened. In the event of a malfunction, the corresponding components should always be replaced in full!

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The following table on page 24 is a guide to rectifying common malfunctions:

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If a function fault is detected during operation, which the user cannot rectify with the use of the trouble-shooting table aid, the specialized staff for maintenance and repair of the appropriate owner is to be called.

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The user may **under no circumstances** try to eliminate defects at the electrical equipment!

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Dangerous voltage, electric shock possible!

Any work on the electrical equipment may only be performed by service personnel, the motor manufacturer or qualified and authorized electricians as required by the applicable VDE rules and safety guidelines!

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Trouble-shooting table aid

Problem	Possible causes	Solution
The patient lift will not go UP or DOWN	<ul style="list-style-type: none"> ■ The emergency off button is depressed ■ The connecting cable is disconnected ■ Battery empty ■ The handset is damaged 	<ul style="list-style-type: none"> ☞ Check, remove if necessary ☞ Check the cable connection ☞ Check the rechargeable battery pack and recharge the battery if necessary ☞ Inspect the keyboard for damage (visual inspection) ☞ Fuse broken - replace fuse ☞ Have keyboard checked by qualified person.
The battery pack will not recharge	<ul style="list-style-type: none"> ■ The fuse is broken. ■ The battery charger is defective 	<ul style="list-style-type: none"> ☞ Check whether the battery charger fuse is OK (Internal 10A fuse) ☞ Check battery charger and replace if necessary
Lift can only be raised or lowered slowly or with difficulty	<ul style="list-style-type: none"> ■ Battery not sufficiently charged 	<ul style="list-style-type: none"> ☞ Check state of battery pack and recharge if necessary

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ACCESSORIES

Lifting accessories, swivel arm attachments in various widths, back harnesses, and seat slings etc., according to the relevant list, are available for all models of the Burmeier patient lift ranges.

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MEASUREMENTS AND DIMENSIONS

Width
Length
Height
Weight base
Weight lift arm
Base: mechanically spreadable using foot pedal
Turning circle (expanded)
Safe working load
Castor set (rear castors locked / dual castors at front)
Noise level during adjustment
Lift actuator (with mechanical emergency lowering)
Power supply (battery system)

Controls:

CBJ H with 1 channel, EMERGENCY OFF and electrical emergency lowering (IPX4)
duty cycle: 10%

Operation:

Hand control with coiled cable

Battery charger: external

All lifts in the range belong to type class B, protection class II.

ELECTRICAL DATA
Control unit

Type
Input voltage
Fuse
Output voltage
Output power
Protection category

Battery charger

Type
Input voltage
Max. current input
Fuse
Output voltage
Output power
Protection category

Lift actuator



Type (alternatively)
Input voltage
Duty cycle (Intermittent duty ON/OFF)
Protection category

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
Item.	Description	Dimension
1	Highest point	–
2	Highest point at maximum range	–
3	Lowest point	–
a	Maximum range at 600 mm height	695
b	Maximum range from base	695
c	Range measured from base with legs spread to 700 mm	220
k	Max. height of swivel bar attachment	1970
l	Min. height of swivel bar attachment	880
m	Lifting range	1090
p	Max. inner width	1050
q	Inner width at maximum range	965
r	Minimum inner width	540
x	Minimum clearance between wall and swivel bar attachment at lowest lift height (base legs spread)	260
y	Minimum clearance between wall and swivel bar attachment at maximum range (base legs spread)	185
z	Minimum clearance between wall and swivel bar attachment at greatest lift height (base legs spread)	610

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TABLES FOR MEDICAL ELECTRICAL EQUIPMENT

General information		
1	Guidelines and Manufacturer's Declaration – Electromagnetic Emissions	
2	The lift is intended for use in the electromagnetic environment described below. The client or user of the lift should ensure that it is used in such an environment.	
3	Spurious emission tests	Compliance
4	HF emissions To CISPR 11	Group 1
5	HF emissions To CISPR 11	Class B
6	Harmonics emissions according to IEC 61000-3-2	Class A
7	Emissions of voltage fluctuations/flicker acc. to IEC 61000-3-3	Not applicable

Interference immunity tests	IEC 60601 - test limit	Upper compliance limit
Electrostatic discharge (ESD) according to IEC 61000-4-2	+/- 6 kV contact discharge +/- 8 kV air discharge	+/- 6 kV contact discharge +/- 8 kV air discharge
Short, transient electrical disturbance / bursts according to IEC 61000-4-4	+/- 2 kV for network cables +/- 1 kV for input and output cables	+/- 2 kV for network cables +/- 1 kV for input and output cables
Surges according to IEC 61000-4-5	+/- 1 kV transversal voltage +/- 2 kV longitudinal voltage	+/- 1 kV transversal voltage +/- 2 kV longitudinal voltage
Voltage dips, short interruptions and fluctuations in the supply voltage according to IEC 61000-4-11	<5% UT (>95% dip in UT) for half a period 40% UT (60% dip in UT) for five periods 70% UT (60% dip in UT) for 25 periods 5% UT (>95% dip in UT) for 5s	<5% UT (>95% dip in UT) for half a period 40% UT (60% dip in UT) for five periods 70% UT (60% dip in UT) for 25 periods 5% UT (>95% dip in UT) for 5s
Supply frequency magnetic field (50/60Hz) according to IEC 61000-4-8	3 A/m	3 A/m

Non life-supporting systems

Interference immunity tests	IEC 60601 - test limit	Upper compliance limit
Conducted HF interference according to IEC 61000-4-6	3V _{eff} 150 kHz to 80 MHz	3V
Radiated HF interferences according to IEC 61000-4-3	3V _{eff} 80 MHz to 2.5 GHz	3 V/m

Protection distances from wireless telecommunications equipment

Nominal output of transmitter W	Protection distance dependent on the transmission frequency		
	150 kHz to 80 MHz $d=P \cdot \exp(0,5 \cdot 3,5/V1)$	80 MHz to 800 MHz $d=P \cdot \exp(0,5 \cdot 3,5/E1)$	800 MHz to 2.5 GHz $d=P \cdot \exp(0,5 \cdot 7/E1)$
0.01	0.12 m	0.12 m	0.24 m
0.1	0.37 m	0.37 m	0.74 m
1	1.17 m	1.17 m	2.34 m
10	3.69 m	3.69 m	7.38 m
100	11.67 m	11.67 m	23.34 m

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APPLIED STANDARDS AND GUIDELINES

Medical Products Directive
Repeat testing ... of medical electrical units
Technical aids for handicapped persons
Handicapped friendly construction
Risk analysis for medical products
Hoists for the transfer of disabled persons
Safety for electro-medical equipment
Electromagnetic tolerances
Classified as class I active medical product (in accordance with MPG § 13).
Type-tested by an independent testing institute

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DISPOSAL INSTRUCTIONS

- ◆ This lift - provided it is electrically adjustable - is classified as a commercially used electrical device (b2b) in accordance with the WEEE Guidelines 2002/96/EC (law governing electrical equipment):
- ◆ Any electrical components of this lift (actuators, control devices, handsets, etc.) that are replaced must be treated as electrical waste in accordance with the WEEE Guidelines and must be properly disposed of.
- ◆ For lifts brought into operation after 13 August 2005 (see the information given on the type plate under "Series"), the operator is legally obliged to send the electrical components directly to the manufacturer and not to dispose of them at municipal waste collection points. BURMEIER and its service and sales partners will take these components back.
- ◆ The return of these components is covered by our general terms and conditions.
- ◆ Individual rechargeable batteries that are no longer of use must be properly disposed of in accordance with battery regulations and do not belong in the household waste.
- ◆ The operator must ensure that all components that are to be disposed of are not infectious or contaminated.
- ◆ If the patient lift is to be scrapped, the synthetic and metallic parts are to be separated and disposed of properly.
- ◆ Gas springs, if present, are under high pressure! They must first be depressurised according to the manufacturer's instructions before disposal. This information can be obtained upon request from the gas spring manufacturer (see type plate).
- ◆ Gas springs and hydraulic units, if present, are made predominantly of metal and plastic. Before disposal, the oil contained in these components must be emptied out in accordance with the manufacturer's instructions and disposed of properly.
- ◆ If you have any queries you can contact your local municipal waste company or our service department.

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EC-DECLARATION OF CONFORMITY

We,
Burmeier GmbH & Co. KG
Pivitsheider Strasse 270
32791 Lage / Lippe
Germany

hereby declare that the product named below:


AGILE Homecare Patient Lift

complies with the regulations of the EC Directive 93/42/EEC for Medical Products.

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No part of this manual may be reproduced with-
out the prior written permission of the publisher!

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All rights reserved.
Subject to technical changes!

Stand: 01.12.2006

Updated: December 01, 2006